

TMEM141-transmembrane protein 141 Gene

Reference: PTXBC007834

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 141		
For research use only.			

GENE DESCRIPTION	transmembrane protein 141
GENE ID	85014
ORIGIN SPECIES	Human
ACCESSIONS	BC007834
VECTOR	pDONR223
NCBI SYMBOL	TMEM141
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgaactgggtctgtcccggtggacgacccgtggctgccaagcaccgggactcgggagtagtccgcatgccagtcacagccttcataagggcggtttcaccttcgtcacaggcaccggcatggccttggcttcagatgttcattcagaggaagtccatacccttgca
gtggagcctcctagtgccgtggttcaggtctgtgtcagctacgggtgacgagagtgagtcggagaaatgaacaacctctggctcttcctggagaccgggcagctcccaagacaggagcacagatcagagaagctag