

C12orf62-chromosome 12 open reading frame 62 Gene

Reference: PTXBC007849

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C12orf62; cytochrome c oxidase assembly protein COX14; COX14 cytochrome c oxidase assembly homolog; cytochrome c oxidase assembly homolog 14; COX14, cytochrome c oxidase assembly factor		

For research use only.

GENE DESCRIPTION	chromosome 12 open reading frame 62
GENE ID	84987
ORIGIN SPECIES	Human
ACCESSIONS	BC007849
VECTOR	pDONR223
NCBI SYMBOL	C12orf62
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccaaatggcaagcagctagctgacattggctataagaccttctctacccatgatgcttctcactgtgtatgggggtacctctgcagtgccgagctaccactattccagtgccgcagggccagcgcagggccagaagaacagaagacctcaggaatcatgtag