

MGC14425-hypothetical protein MGC14425 Gene

Reference: PTXBC007866

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein MGC14425
GENE ID	84989
ORIGIN SPECIES	Human
ACCESSIONS	BC007866
VECTOR	pDONR223
NCBI SYMBOL	MGC14425
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacccccgtcgcagcttcgagccggtccgctctccagcgtcgaacgtgcctcgtccgaccgccacacaggaaccgcttaccacacagctctccccgctctctaccgcatagctgtcgtcgcgaagcggccgtgcctcctccagtcgaggggaaccgatga
aacctcactcttaccggcgtctatgctgagggagcggaccgggacacagcagcggacccgaaagagcgcagactcgggacgaaccggcgtctgccccggacacagcagctcgggcctccccgaaaacactccttggactcccagattcgcagccttgtgtgcagcgc
acacaagaaaactga