

ZNF625-zinc finger protein 625 Gene

Reference: PTXBC007868

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein 625		
For research use only.			

GENE DESCRIPTION	zinc finger protein 625
GENE ID	90589
ORIGIN SPECIES	Human
ACCESSIONS	BC007868
VECTOR	pDONR223
NCBI SYMBOL	ZNF625
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggagagagactcttagaaagtaagaagatcatcagcatggagaaatttgaccaggtccagatgacatgctgaagaagaaacccccagtaaaatcatgtggagaagtcagcgtgggtcatgcatcccttaataggcaccagagctgacactggacacaagccatagat
atcaggaaatggacagaagccataataatgtacatactgtaagaagccttcagtgatctccctacttccgaacacatgaatgggtcacactgggggaaaccttatgattgtgaggaaatgggaaagccttatttccgltcaagcattcgaagacacaggataatgcacagtgagat
ggaccctacaatgtaacttttgggaaagccttgatggtcagtttgcacacacgaactcacactggagagaaccatgatgaatgaacagtggtgaaagcctttagtcattctgtagccttcgaatacatgaagaactcacactggagagaagccttatgaatgcagtg
agtgtgggaaagcattccatagttccatgcctccatgcacataaaataactcacactgggagaagccgtatgaatgaacagtggtgaaagccttgggttcttaattcgttcgatatcatgaagaactcacactggagagaagccttatgaatgaacagtggtgaaagccttc
agatctgcctcgacactcgaacacatggaaggactcacactggagagaagccttatgaatgaacagtggtgaaagccttgggtgctcgcagcgttaaaatccatgaagactcacactggagaaagcctttagctccaacactcgaaggccaaggcagagaagattgctt
aa