



Reference: PTXBC007897

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	P97; RPN1; TRAP2; 26S proteasome non-ATPase regulatory subunit 2; 55.11 protein; TNFR-associated protein 2; proteasome (prosome, macropain) 26S subunit, non-ATPase, 2; protein 55.11; tumor necrosis factor type 1 receptor-associated protein 2; proteasome 2		

For research use only.

GENE DESCRIPTION	proteasome (prosome, macropain) 26S subunit, non-ATPase, 2
GENE ID	5708
ORIGIN SPECIES	Human
ACCESSIONS	BC007897
VECTOR	pDONR223
NCBI SYMBOL	PSMD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]