

SPIB-Spi-B transcription factor (Spi-1/PU.1 related) Gene

Reference: PTXBC007921

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SPI-B; transcription factor Spi-B; Spi-B transcription factor (Spi-1/PU.1 related); Spi-B transcription factor		

For research use only.

GENE DESCRIPTION	Spi-B transcription factor (Spi-1/PU.1 related)
GENE ID	6689
ORIGIN SPECIES	Human
ACCESSIONS	BC007921
VECTOR	pDONR223
NCBI SYMBOL	SPIB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgcccctggaggctgacagctcagcggccacactcagctgtctgtaccagatggcgtcttctatgacctggacagctgcaagcattccagctaccctgattcagagggggctcctgactccctgtgggactggactgtggccccctgtccagccacccctatgaagcctt
cgacccggcagcagccgcttttagccaccccccagctgcccagctctgtctacgaacccccacactacagccctgagggaaacctggaactggccccagcctggaggccccggggcctggcctccccgataccccagggagaacttcgtagccagacccctgttccccggcat
atgccccgtacccagccctgtgtatcagaggaggaagactaccgttgacagccctgcccctggaggctcgcgacagcagtgatgggacccctcgtggctggccccgaggggaaggatccgaggcaggactcgcaagaagctgcctgtaccagttcctgtggggcta
ctgacgcgcgggacatgctgagtgctgtgtgggtggagccagcgcggcgtcttcagttctctccaagcacaaggaactcctggcgccgctggggccagcagaaggggaaccgcaagcgcgtacgtacacagaagctggcgcgcgcctccgaaactacgccaag
accggcgagatccgaaggcaagcgaagctcacctaccagttcagacagcgcgtgctgctgacgtccgcccggcctga