

C21orf119-chromosome 21 open reading frame 119 Gene

Reference: PTXBC007928

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	chromosome 21 open reading frame 119		

For research use only.

GENE DESCRIPTION	chromosome 21 open reading frame 119
GENE ID	84996
ORIGIN SPECIES	Human
ACCESSIONS	BC007928
VECTOR	pDONR223
NCBI SYMBOL	C21orf119
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggcagcgggacactccacgacacccccccccggctgcaggggtttgggtccacggggcgcatcttgatccagtggtcttgggtgctcttggccggaagaactccagaccttcactcctggcctggccgatccacgaagactcgggccccggcggtccaagcag
cacaatga