

GAP43-growth associated protein 43 Gene

Reference: PTXBC007936

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	nerve growth-related peptide GAP43; B-50; PP46; neuromodulin; axonal membrane protein GAP-43; calmodulin-binding protein P-57; neural phosphoprotein B-50; neuron growth-associated protein 43; protein F1; growth associated protein 43		

For research use only.

GENE DESCRIPTION	growth associated protein 43
GENE ID	2596
ORIGIN SPECIES	Human
ACCESSIONS	BC007936
VECTOR	pDONR223
NCBI SYMBOL	GAP43
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgtgtgtatgagaagaacacaggttgaaaaaatgatgacgacaaaagattgaacaagatggtatcaaacagaagataaagctcataaggccgcaacaaattcaggctagcttccgtggacacataacaaggaaaaagctcaaaggagagaagaagatgatgtcca
agctgtgaggtgagctaaataagaagatgaagccctgttgcgatgggtggagaagaaggagagaaggcaccactactgccgaagcagcccgccactggctccaagcctgatgagcccgcaagcaggagaaactctccgaggagaagaaggagggtgatgc
tgccacagagcaggcagccccaggtctctgcatctcagaggagaaggccgctcagctgagacagaagtgccactaaagcttcactgataactgccgtctcccaaggctgaagatgcccgcccaaggagagagcctaacaagccgatgtcctgtctgtcactgtctgtc
gctgccaccacccctgccgagagatgctgtgccaaggcaacagccagcctccaacggagactggggagagcagccaagctgaagagaacatagaagctgtatgaaaccaaactaaggaaagtgcccgagagagggtaagaagagggaacctgaggctgacca
agaacatgcctga