



Reference: PTXBC007981

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	serine/threonine-protein kinase STK11; PJS; hLKB1; liver kinase B1; polarization-related protein LKB1; renal carcinoma antigen NY-REN-19; serine/threonine-protein kinase 11; serine/threonine-protein kinase LKB1; serine/threonine kinase 11		

For research use only.

GENE DESCRIPTION	serine/threonine kinase 11
GENE ID	6794
ORIGIN SPECIES	Human
ACCESSIONS	BC007981
VECTOR	pDONR223
NCBI SYMBOL	STK11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]