

## ARHGEF18-rho/rac guanine nucleotide exchange factor (GEF) 18 Gene

Reference: PTXBC008016

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | P114-RhoGEF; rho guanine nucleotide exchange factor 18; 114 kDa Rho-specific guanine nucleotide exchange factor; Rho-specific guanine nucleotide exchange factor p114; Rho/Rac guanine nucleotide exchange factor (GEF) 18; SA-RhoGEF; p114RhoGEF; septin-assoc |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | rho/rac guanine nucleotide exchange factor (GEF) 18                                               |
| <b>GENE ID</b>          | 23370                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC008016                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | ARHGEF18                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgtccaggcagatgcagagatgcacctagacgttgcagcaagtggacaagtggccgctgtgcgggcccctcgctgtagtgcagctgtgcagcttacggcttcctcctggagggtggaggaggaggtgtgggcagcatcaaaagtgctgggacatccagggtggtgaga  
tcacccacgatccagctccggtggagaaaggcccatgtcaagcctgttctgcacccaagcaatgggtgtaggactgggtcctgctga