



Reference: PTXBC008025

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C16orf21; JFP7; WD repeat-containing protein 24; WD repeat domain 24		
For research use only.			

GENE DESCRIPTION	WD repeat domain 24
GENE ID	84219
ORIGIN SPECIES	Human
ACCESSIONS	BC008025
VECTOR	pDONR223
NCBI SYMBOL	WDR24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr