

LOC158435-hypothetical protein BC008050 Gene

Reference: PTXBC008050

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC158435; uncharacterized protein LOC158435		

For research use only.

GENE DESCRIPTION	hypothetical protein BC008050
GENE ID	158435
ORIGIN SPECIES	Human
ACCESSIONS	BC008050
VECTOR	pDONR223
NCBI SYMBOL	LOC158435
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatttcagcaggctggtggaggcagagcagagaaggcgggagaggccagcgaaggcctgggaagggtggacccacaactcgtctctaataagtaaccagcaatgcgagggcacacagggagacccctagaggctgccccagggcaggagtcacctgtgaactg
a