

## RP5-1077B9.4-invasion inhibitory protein 45 Gene

Reference: PTXBC008068

|                |                                                                                                                                                     |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                 | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | IIP45; migration and invasion-inhibitory protein; IGFBP2-binding protein; invasion-inhibitory protein 45; migration and invasion inhibitory protein |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | invasion inhibitory protein 45                                                                    |
| <b>GENE ID</b>          | 60672                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC008068                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | RP5-1077B9.4                                                                                      |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggtggaggctgaggactggcacagctgcggctgctcaatctggagctctgaggcagctgtgggtggggcaggatgctgtgcggcggtcagtgccaggcgagcctcggagtcaggcctggaatccagcagcagctacaactcagagactccatcgacccagagagctcctc  
aacttccttgagcactcctgcccacggggccggctcctcctcctgctggggccaccagatgctgtgaggggacccctcctgctgagtgccagatcgggggtggcctctctccacctgccaaftgccagcaccaggagtccttggccgaccgagacccactcagcaccctcctgggc  
acctcaagcctgaggggaccagagccctcaggaggctgggtgatccagaccagaggagacacagctcaggtccatcctgctcaacagagcaagctgtccaagccagggtgaccttctctgaggagctcaggtctcagagagagctggcgctcaggccatcctggg  
ctatgactggattgcagggtctctggacaccagctctccatcaccagccagcctgaggcctctctcctcaagctgcaggaggttgcggaaaccaacaaggaggagtgatctgcagccatcctgaacccagttgccaggcctgcgtgagagcagtgggcagcggcgtggaggaagacc  
atgaatgcgtgactgttaccgtgtcaaccggcgcctgttcccggtgctgtggtatcccggtaccocctgcgcctgtgcaggacacccgcgagaccagcaggggccctggggaccctggcgagccagcgcagctcagggtgagcatcccctgtgatcctggagccccgcaccggt  
accacatccacggcgaaagagctttgacgccttgacacactggccctgccccggcactgctgctgggctgggacatttttctccgaagtctgagaaaagctcagccccaggaaactggagctctggtcctctgtctcctgaggccagcaccagaagctgtccggcaccagca  
gcccttttaccggcctcaccatgcatgctgccccgaccccgacccctgctcagtgccccaggctccctggccccacgtccacggcagaagccctga