



Reference: PTXBC008077

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SRPR; Sralpha; signal recognition particle receptor subunit alpha; DP-alpha; SR-alpha; SRP-alpha; docking protein alpha; signal recognition particle receptor (docking protein); SRP receptor alpha subunit		

For research use only.

GENE DESCRIPTION	signal recognition particle receptor (docking protein)
GENE ID	6734
ORIGIN SPECIES	Human
ACCESSIONS	BC008077
VECTOR	pDONR223
NCBI SYMBOL	SRPR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]