

## CYTH3-cytohesin 3 Gene

Reference: PTXBC008191

|                |                                                                                                                                                                                                                                        |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                    | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                   | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | ARNO3; GRP1; PSCD3; cytohesin-3; ARF nucleotide-binding site opener 3; PH, SEC7 and coiled-coil domain-containing protein 3; general receptor of phosphoinositides 1; pleckstrin homology, Sec7 and coiled-coil domains 3; cytohesin 3 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | cytohesin 3                                                                                       |
| <b>GENE ID</b>          | 9265                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC008191                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | CYTH3                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgaaccgcggcatcaacgagggcggggacccctcctgagagctgctgaggaaattgtatgagagcattaagaacgagccatttaagatcccgaggagcagcgggaacgacctgacccacaccttctcaaccccgaccgcgagggtggtcctgaagctggaggggcgtgtgaa  
gacctggaagcgcgggtgttcatctgaccgataactgcctctattacttgaatacacaagataaggagccagggaatcccggttgaaaacctcagcatcaggaggtggaggaccccgaaacccaactgtttgagctctacaatccagccacaaggcagggtcat  
caaggcctgtgaagacggaggccgacggccgctgttagaggggaacctgtgtgtaccggaatcagccccgagccggaggagaaggaggagtgatgaaatccatcaaaagccagtcagcagagatccctctatgacatgttggaacgaggaaacgaaggattgccaata  
aaaaatag