

HEBP2-heme binding protein 2 Gene

Reference: PTXBC008205

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6ORF34B; C6orf34; PP23; heme-binding protein 2; placental protein 23; heme binding protein 2		

For research use only.

GENE DESCRIPTION	heme binding protein 2
GENE ID	23593
ORIGIN SPECIES	Human
ACCESSIONS	BC008205
VECTOR	pDONR223
NCBI SYMBOL	HEBP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgagccgctccagccagacccggggcgccgagggcgccgccaagctgtggagacgcccggctggaagggccggaggagcggccggccggaggcgggaagttagatccgacactatggaccgcaagtgggtcagcagctcgtggagctatggact
gggattcagccatccagacgggctttacgaaactgaacagctacatcaaggcgggagagagatgaaataaagatgacagctccagtgacaagctacgtggagccgtgggtcaggtccttttagtgagctaccattaccattccctgtatattccctctgaacagcaattgatcca
cccaggcccttagagtcagatgtctcattgaagatagagccgaaatgactgtgtttgtacggctcttcgatggatttttagtgcccaaaagaatcaagaacaacttttgacattagcaagcattttaagggaagatggaaaagtgttcgatgagaaggtttactacactgcaggctacaacagtc
tgtcaaatgtcttaataagaataatgaagtgtgttgattcaaaaaatgaaccacaaagaacgaatga