

C9orf46-chromosome 9 open reading frame 46 Gene

Reference: PTXBC008212

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein C9orf46; C9orf46; AD025; MDS030; PLG-RKT; Plg-R(KT); plasminogen receptor (KT); 5033414D02Rik; plasminogen receptor, C-terminal lysine transmembrane protein; plasminogen receptor with a C-terminal lysine		

For research use only.

GENE DESCRIPTION	chromosome 9 open reading frame 46
GENE ID	55848
ORIGIN SPECIES	Human
ACCESSIONS	BC008212
VECTOR	pDONR223
NCBI SYMBOL	C9orf46
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtttatatttcaaatctatgaatgaaagcatgaaatcaaaaggagttcattgaatgctcgaactcagctggaaggcagctcatcagagtgaaatgagggaagacaatggccatgcagattgcgtggtctcgggaattcctcaaatatttggaaactttttggcct
tgcagccatctttaaagcgtggagcgattaaaaaaagaagccagcttctggtcccgattgttcattaaagctttatcctcacctaccagtatgacttggcctatggaaccctttagaagaatgaaagggtgaagctgaggacatactggaacagaaaagagtaaatgcagctgcaa
gaggaatgatcactttgaaagcattgaaaaagccagaaaggacagagtagattctcatagacaaatga