

ABHD11-abhydrolase domain containing 11 Gene

Reference: PTXBC008251

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein ABHD11; PP1226; WBSCR21; Williams Beuren syndrome chromosome region 21; abhydrolase domain-containing protein 11; alpha/beta hydrolase domain-containing protein 11; abhydrolase domain containing 11		

For research use only.

GENE DESCRIPTION	abhydrolase domain containing 11
GENE ID	83451
ORIGIN SPECIES	Human
ACCESSIONS	BC008251
VECTOR	pDONR223
NCBI SYMBOL	ABHD11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggccatcaacatgcagatgagctgccccgtcccggtgccccgaaactggcgatgaacagctcagttctgtcatccaggacatggccgtgcccagcacctgtcactaacctggtagaggtagacggcgcttcgtgtgagggtgaactggatgccctgaccagcacct
agacaagatcttgctttcccacagaggcaggagtcctaccctgggccaacactctttctcttgggtgaaactccagttcgtgcatccagccaccacccctgagattatgcggctcttcctcgggccagatgcagacgggtgccgaacgctggccactggatccacgctgaccgcca
caggacttcatagtgccatccgaggtctctgtctaa