

PTGES-prostaglandin E synthase Gene

Reference: PTXBC008280

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MGST-IV; MGST1-L1; MGST1L1; MPGES; PGES; PIG12; PP102; PP1294; TP53I12; mPGES-1; MGST1-like 1; glutathione S-transferase 1-like 1; microsomal glutathione S-transferase 1-like 1; microsomal prostaglandin E synthase-1; p53-induced apoptosis protein 12; p53-		

For research use only.

GENE DESCRIPTION	prostaglandin E synthase
GENE ID	9536
ORIGIN SPECIES	Human
ACCESSIONS	BC008280
VECTOR	pDONR223
NCBI SYMBOL	PTGES
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgccacagcctggtgatgagcagcccgccctccggcctctctgctctgcagcagcgtgctgtgcatcaagatgtacgtggtggccatcatcacgggccaagtggagctgcggagaaggcctttgccaaaccccgagatgccctgagacacggagggccccagttattgca
ggagcgaccccgacgtggaacgctgctcagggccaccgggaacgacatggagaccatctaccctctctttctgggcttcgtctactcctttctgggtcctaaccctttgtcgctggatgcacttctgtcttctcgtggcgctgtggcacacacgtggcctacctggggaagct
gcgggcacccatccgtccgtgacctacacctggccagctccctcgcgctccatggcctctgcagatcctctgggaagcgcccgccacctgtga