



Reference: PTXBC008281

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GMPR 1; GMPR1; GMP reductase 1; guanine monophosphate reductase; guanosine 5'-monophosphate oxidoreductase 1; guanosine monophosphate reductase 1; guanosine monophosphate reductase		

For research use only.

GENE DESCRIPTION	guanosine monophosphate reductase
GENE ID	2766
ORIGIN SPECIES	Human
ACCESSIONS	BC008281
VECTOR	pDONR223
NCBI SYMBOL	GMPR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]