

PCNXL2-pecanex-like 2 (Drosophila) Gene

Reference: PTXBC008300

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PCNXL2; pecanex-like protein 2; pecanex homolog 2 (Drosophila)		
For research use only.			

GENE DESCRIPTION	pecanex-like 2 (Drosophila)
GENE ID	80003
ORIGIN SPECIES	Human
ACCESSIONS	BC008300
VECTOR	pDONR223
NCBI SYMBOL	PCNXL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgtttttgtggttctgtgtctgtggaataacctcggtgtttacagtgtggccggagcgtcttggctgccgccctgctccacgcagtgtcttcagtcagtgagggaaccgtggagcatgcaacacatcccgcaactgtttcggccttctgtggccttctgtgcccttttaccatctgagccgtcagagcagtgacccatctgtactcatgtcttccaatgcaggctgtttcctaaattttacatcaaaatctggcagagtcagctgtgacctctcccaagaagatgaaagattcagtgacggatgtcttaaaatgggatctcatgtctgcgcagtggttgcctctcatttgagtcagcgcagcagcagctgtattcctgtcattgcgacattctcagcatcgtgtgttgccttggctggagccgtgggtgtttgtaacacattacgtgctccctcagctccgcaagcatcatccctggatgtggatttcacaccccatctcaaaaacaaagagtatcatcaacgggaagtgagagatgttgcctttaaagtggttcgaaagactctatgttggctcagtggtttgaaaaatacatcttgaaccagcgctaatttgaatgccctcactattgatgcatgttttaataagcaatcaccggagacttgggtaccacatgttttctgatgatcattggcatgaagctgttgcggacatcattctgcaacccgtttaccagtttataacttgagcttcactgtcatcttttccacttgaactacaaagatatcagagagcttctactggatttctcatgtgtccattttatagcaagggttttctgggtacgtgagggcaaaaataa