



Reference: PTXBC008368

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C12orf11; GCT1; INTS13; Mat89Bb; NET48; SPATA30; protein asunder homolog; asunder, spermatogenesis regulator homolog (Drosophila); cell cycle regulator Mat89Bb homolog; germ cell tumor 1; integrator complex subunit 13; sarcoma antigen NY-SAR-95; spermatoge		

For research use only.

GENE DESCRIPTION	chromosome 12 open reading frame 11
GENE ID	55726
ORIGIN SPECIES	Human
ACCESSIONS	BC008368
VECTOR	pDONR223
NCBI SYMBOL	C12orf11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr