

CHP-calcium binding protein P22 Gene

Reference: PTXBC008373

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	calcium-binding protein CHP; CHP; SLC9A1BP; Sid470p; p22; p24; calcineurin B homologous protein 1; EF-hand calcium-binding domain-containing protein p22; SLC9A1 binding protein; calcineurin B homolog; calcineurin B-like protein; calcineurin homologous pro		

For research use only.

GENE DESCRIPTION	calcium binding protein P22
GENE ID	11261
ORIGIN SPECIES	Human
ACCESSIONS	BC008373
VECTOR	pDONR223
NCBI SYMBOL	CHP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagtctcactctgtcaccaggctggagtcagtgccgtgatcttgctcactgcaacctctgcctctgggttcaagcaattctcccactcagcctcccaagtagctgggattacagacgtgtgccaccatacctggtaattttgcattttttagtgagaggaggagttcaccatgttgccaggttggtcttga