



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CXCL5-chemokine (C-X-C motif) ligand 5 Gene

Reference: PTXBC008376

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ENA-78; SCYB5; C-X-C motif chemokine 5; chemokine (C-X-C motif) ligand 5; epithelial-derived neutrophil-activating protein 78; neutrophil-activating peptide ENA-78; neutrophil-activating protein 78; small inducible cytokine subfamily B (Cys-X-Cys), member		

For research use only.

GENE DESCRIPTION	chemokine (C-X-C motif) ligand 5
GENE ID	6374
ORIGIN SPECIES	Human
ACCESSIONS	BC008376
VECTOR	pDONR223
NCBI SYMBOL	CXCL5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]