

EPB41L3-erythrocyte membrane protein band 4.1-like 3 Gene

Reference: PTXBC008377

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	4.1B; DAL-1; DAL1; band 4.1-like protein 3; differentially expressed in adenocarcinoma of the lung protein 1; erythrocyte membrane protein band 4.1 like 3		

For research use only.

GENE DESCRIPTION	erythrocyte membrane protein band 4.1-like 3
GENE ID	23136
ORIGIN SPECIES	Human
ACCESSIONS	BC008377
VECTOR	pDONR223
NCBI SYMBOL	EPB41L3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagtcgcaagtgatactctcgatggatcagaatatactgtgatgtagagaacgctccagagacaagtgctgttgataaagtgtgtgaacactgaacttgctagagaaagactactttgggcttacgtatcgagatgctgaaaaccagaagaattggttgaccctgctaaggaaat
aaaaaacagggtcgaagtgtgcttgccactttcatttaattgtgaaattttaccaccagacctgcccactatctgaagatataccagggtactactctgcttcagttgcgagatgacatcgtgtccggaaggctgccctgctcctttgtaccctggccttgcctggctcctaa