

IMPA1-inositol(myo)-1(or 4)-monophosphatase 1 Gene

Reference: PTXBC008381

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IMP; IMPA; inositol monophosphatase 1; D-galactose 1-phosphate phosphatase; IMP 1; IMPase 1; inositol(myo)-1(or 4)-monophosphatase 1; inositol-1(or 4)-monophosphatase 1; lithium-sensitive myo-inositol monophosphatase A1; myo-inositol monophosphatase 1; te		

For research use only.

GENE DESCRIPTION	inositol(myo)-1(or 4)-monophosphatase 1
GENE ID	3612
ORIGIN SPECIES	Human
ACCESSIONS	BC008381
VECTOR	pDONR223
NCBI SYMBOL	IMPA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgaccttggcagggaatgcattatgcagtaactctagcaagacaagctggagaggtgtgtgtgaagctataaaaaatgaatgaatgttatgctgaaaagtctccagttgatttgtaactgctacggacacaaaagtgaataatgcttatctctccataaaggaaaagtatc
catctcacagtttcaatgggtgaagaatctgtggcagctggggaataatgtatcttaaccgacaacccacatggatcattgacccattgatggaacaactaacttgcacatagatttctttgtagctgtttcaatggccttgcgtgtaataaaaaagatagaattggagttgtgacagttgtgtg
gaaggcaagatgtacactgccagaaaaggaaaagggtgctttgtatgggtcaaaaactacaagctctcacaacaagaagatattaccaatctctctgtgactgagttggcctcttcagaacaccagagactgtgagaatggttcttctaataatggaaggctttttgcatcctgttcattgg
gatccggagttgtggaacagcagctgttaatatgtgcctgtggcaactggcggagcagatgcattattgaatgggaattcactgctgggagttgcaggagctggcattattgtactgaagctggtggcgtgctaattggtacaggtggaccattgattgatgtcacgaagagta
tgctgcaataatagaatattagcagaaaggatagctaaagaattcagggttatacctttgcaacgagacgacgaagattaa