

ELA3A-elasticase 3A, pancreatic Gene

Reference: PTXBC008383

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ELA3A; ELA3; chymotrypsin-like elastase family member 3A; elastase 1; elastase 3A, pancreatic; elastase IIIA; elastase-3A; protease E; chymotrypsin like elastase family member 3A		

For research use only.

GENE DESCRIPTION	elastase 3A, pancreatic
GENE ID	10136
ORIGIN SPECIES	Human
ACCESSIONS	BC008383
VECTOR	pDONR223
NCBI SYMBOL	ELA3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgccagctgctcagttccctctctctgtggcgttgcctcaggtatggccaccttctctcactcttccagccgctgttccatggtgaggatgcggtccctacagctggccctggcaggttccctgcagtatgagaaaagtgaagcttctaccacacgtgtgcccgtagcct
catcgcctccgaltgggttgactgccggccactgcatctgagggatctgacctaccaggtgtgtgtgggtgagtacaaccttgcctggaaggaggcccccagcaggtgatcccatcaactctgaggagctgtttgtgcatccactctggaaccgctcgtgtgtggcctgtggcaatg
acatgccctcatcaagctctcagcagcggccagctgggagatgccgtccagctcgcctcactccctccgctggtgacatctctcccaacaagacacctgctacatcaccgctggggcgtctctataccaatgggccaactcccagacaagctgcagcagcccgctgccgtg
gtggactataagcactgctcaggtggaaactgtgggttccaccgtgaagaaacatggtgtgtgctgaggggtacatccgctccggtgcaacggtgactctggaggaccctcaactgccccacagaggatggtggctggcaggtccacggtgtgaccagctttttctgcttgg
gctgcaacttcatctggaagccacggtgttccactcagctccgcttcatcagctggtgaggagacatagcaagccactag