

ZNF273-zinc finger protein 273 Gene

Reference: PTXBC008386

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HZF9; zinc finger protein 273; zinc finger protein 9; zinc finger protein HZF9		

For research use only.

GENE DESCRIPTION	zinc finger protein 273
GENE ID	10793
ORIGIN SPECIES	Human
ACCESSIONS	BC008386
VECTOR	pDONR223
NCBI SYMBOL	ZNF273
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaccactgacatttagggatgtggccataaaattctctgaggagtggaatgcctggacacttcacagcagaattgtataggaatgtgatgttagataactacagaaacctggctctcctgggtattgctgtcttaagccagacctgatcactgtctggagcaaggaaaagagcc
ctgcaatatgaagagacatgcgatggtagccaaacccaggttggagtgacgtggcgtga