



Reference: PTXBC008406

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD36 molecule; leukocyte differentiation antigen CD36; CD36 molecule (thrombospondin receptor); CD36 antigen (collagen type I receptor, thrombospondin receptor); BDPLT10; CHDS7; FAT; GP3B; GP4; GPIV; PASIV; SCARB3; platelet glycoprotein 4; GPIIIB; PAS IV;		

For research use only.

GENE DESCRIPTION	CD36 molecule (thrombospondin receptor)
GENE ID	948
ORIGIN SPECIES	Human
ACCESSIONS	BC008406
VECTOR	pDONR223
NCBI SYMBOL	CD36
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]