

TPM3-tropomyosin 3 Gene

Reference: PTXBC008407

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CAPM1; CFTD; HEL-189; HEL-S-82p; NEM1; OK/SW-cl.5; TM-5; TM3; TM30; TM30nm; TM5; TPMsk3; TRK; hscp30; tropomyosin alpha-3 chain; alpha-tropomyosin, slow skeletal; cytoskeletal tropomyosin TM30; epididymis luminal protein 189; epididymis secretory sperm bi		

For research use only.

GENE DESCRIPTION	tropomyosin 3
GENE ID	7170
ORIGIN SPECIES	Human
ACCESSIONS	BC008407
VECTOR	pDONR223
NCBI SYMBOL	TPM3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatggagccatcaagaaaagatgcagatgctgaagttagacaaggagaatgctggtatggcagagcaagctgaagctgagcagaagcaggcagaagaagaagtaaacagctggaggatgagctggcagccatgcagaagaagctgaaaggacagaggatgagctg
gacaagtattctgaagcttgaaggatgccaggaagctggaactggcagagaagaaggctgctgctgctgagctgaggtggcctcttgaaccgttaggatccagctgggttgaagaagagctggacgtgctcaggagcgctggccactgccctgcaaaagctggaagaagctg
aaaaagctgctgagagtgagaggtatgaaggcttattgaaacccggccttaaaagatgaagaaaagatggaactccaggaaatccaactcaaaagctgaagcaccattgcagaagagcagataggaaagtatgaagaggtgctcgttaagtggatcattgaaggagacttg
gaacgcacagaggaacgagctgagctggcagagctaaagtgttctgagctggaggagagctgaagaatgtcaccaacaacctcaagctctttagggctcaggcggagaagtactctcaaaaagaagataaatatgaggaagaatcaagattctactgataaactcaaggaggcaga
gacccgtgctgagtttctgagagatcggtagccaagctggaaaagacaattgatgacctggaagatgagctctatgcccaagaactgaagtacaaggccattagcaggagctggaccacgcccctcaatgacatgacctctatataa