

C15orf15-chromosome 15 open reading frame 15 Gene

Reference: PTXBC008409

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C15orf15; HRP-L30-iso; L30; RLP24; RPL24; RPL24L; TVAS3; 60S ribosomal protein L30 isolog; homolog of yeast ribosomal like protein 24; my024 protein; ribosomal L24 domain-containing protein 1; ribosomal L24 domain containing 1		

For research use only.

GENE DESCRIPTION	chromosome 15 open reading frame 15
GENE ID	51187
ORIGIN SPECIES	Human
ACCESSIONS	BC008409
VECTOR	pDONR223
NCBI SYMBOL	C15orf15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgtatcgaaaagtgttattctgttcggggcccatctatctggacacggcatgatgttcgtccgcaacgattgcaaggtgttcagattttgcaaatctaaatgtcataaaactttaaaagaagcgcaatcctcgcaagtttaggtggaccaaaagcattcggaaagcagctgtaaaaga
gcttacagtggataatcatttgaatttgaaaaacgtagaatgaacctatcaaataccagcgagagctatggaataaaactattgatgcgatgaagagattgaagaaatcaaacagaaagcgccaagctaaatttataatgaacagattgaagaaaaataaagagctacagaaagtcaggat
atcaaaagaagtcagcaaaacatccatcttatccgagccctcttcaggcaaaagggaacagttggaagagaaaaatggtacagcagttacaagaggatgtggacatggaagatgctccttaa