

SRP9-signal recognition particle 9kDa Gene

Reference: PTXBC008443

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ALURBP; signal recognition particle 9 kDa protein; signal recognition particle 9kD; signal recognition particle 9kDa; signal recognition particle 9		

For research use only.

GENE DESCRIPTION	signal recognition particle 9kDa
GENE ID	6726
ORIGIN SPECIES	Human
ACCESSIONS	BC008443
VECTOR	pDONR223
NCBI SYMBOL	SRP9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgcagtagcagacacgtgggagagttcagccgcgtgccgagaagctttacctcgtgaccctatgaaggcagctgtgttctcaatatagcattctgatgggaactgtgtgttaaagtaacagatgatttagttttgtgtataaaacagaccaagctcaagatgtaagaagat
tgagaaattccacagtcactaatgcgacttatggtagccaagggaagcccgcaatgttaccatggaaactgagtga