

SNX3-sorting nexin 3 Gene

Reference: PTXBC008444

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Grd19; MCOPS8; SDP3; sorting nexin-3; sorting nexin 3A; sorting nexin 3		
For research use only.			

GENE DESCRIPTION	sorting nexin 3
GENE ID	8724
ORIGIN SPECIES	Human
ACCESSIONS	BC008444
VECTOR	pDONR223
NCBI SYMBOL	SNX3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggagaccgtggctgacacccggcggctgacaccaaagccgagaaactgaatgacgcctacggacccccagcaactctcagatcagtgtagcaacccgcaaacgggtgggggtcggccggggccgctcaccacttacgaaatcagggtcaagacaaatctctattt
tcaagctgaaagaatctactgttagaagaatacagtgactttgaatggctgcgaagtgaattgaaagagagagcaaggtcgtagttccccgctccctgggaaagcgtttttgcgtcagcttcttttagaggagatgatggaatattgatgacaattttatgaggaaagaaacaagg
ctggagcagtttataacaaggtcgtgctgcatcctctggcacagaacgaacgtgtgtctcacatgtttttacaagatgaataatagataaaagctatactccatctaaataagacatgcctga