

**MLLT11-myeloid/lymphoid or mixed-lineage leukemia (trithorax homolog, Drosophila);
translocated to, 11 Gene**

Reference: PTXBC008445

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AF1Q; protein AF1q; ALL1 fused gene from chromosome 1q; myeloid/lymphoid or mixed-lineage leukemia (trithorax homolog, Drosophila); translocated to, 11; myeloid/lymphoid or mixed-lineage leukemia; translocated to, 11		

For research use only.

GENE DESCRIPTION	myeloid/lymphoid or mixed-lineage leukemia (trithorax homolog, Drosophila); translocated to, 11
GENE ID	10962
ORIGIN SPECIES	Human
ACCESSIONS	BC008445
VECTOR	pDONR223
NCBI SYMBOL	MLLT11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggaccctgtgagtagccagtagcttcttcttctgagatgccatcccagaaactggatctgtcggagctggaaggcctgggtctgtcagatagccacactacaagggtcaagacagcagcgttggcaaatgatgggcaagcaactgcagcagaccaggagaaaaac
cctgaagggtgatggcctcttgagtagcagcacttcaacttctggagagctccattgccagcatccactcttgaactggactgtcttaa