

USE1-unconventional SNARE in the ER 1 homolog (S. cerevisiae) Gene

Reference: PTXBC008455

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	USE1-like protein; vesicle transport protein USE1; D12; MDS032; P31; SLT1; Q-SNARE; SNARE-like tail-anchored protein 1 homolog; protein p31; unconventional SNARE in the ER 1 homolog; unconventional SNARE in the ER 1		

For research use only.

GENE DESCRIPTION	unconventional SNARE in the ER 1 homolog (S. cerevisiae)
GENE ID	55850
ORIGIN SPECIES	Human
ACCESSIONS	BC008455
VECTOR	pDONR223
NCBI SYMBOL	USE1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgtcaggctggagctaaacctgggtgcgctctatcccgtgcgagggcgatggcagcggagaaacgggacccggacgagtggtgcgctggagaaagtacgtggagccctagaggacatgttcaggccctgaaggtccacgcgagcaaacggcctctgaggtgac
aatgaatattctggaagtggtatctgaaggggatctgcaagccgagaagctgacctctctcagagaaagcactggccaaccagtctcgtggccctggcgtgtgccaaccacagccagagagcagtgccccccacaagacgggtgcatctgcagtcacggcggtgaca
ccagcgagatgcggagtgagctactagggcacggactctgcagagcctgagatggacgttaaggagagaactggagtggtggcaggtgccagccagtgagtgagaagcagtcggcagctgagctagacctgctctgcagcgacatcagaacctccaggaagctggcgaagagat
gctaggactggccggagcctcaagaccaataacctggccgccagagtgctcatcaagaaggacaaccagacctgtcacactcactgaaatggcgaccagaacctggagaaactgaagacggagtcagagcgtctggagcagcacgcagaaagtcagtcactggctgctct
gggcatgctcattatgctgtctcatcttattagcatgatcctcttcattcgaatcatgcctaaactaaataa