

FMC1-formation of mitochondrial complexes 1 homolog (S. cerevisiae) Gene

Reference: PTXBC008467

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C7orf55; HSPC268; UPF0562 protein C7orf55; formation of mitochondrial complexes 1 homolog; formation of mitochondrial complex V assembly factor 1 homolog		

For research use only.

GENE DESCRIPTION	formation of mitochondrial complexes 1 homolog (S. cerevisiae)
GENE ID	154791
ORIGIN SPECIES	Human
ACCESSIONS	BC008467
VECTOR	pDONR223
NCBI SYMBOL	FMC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcccttaggtccccggcgacacttttcgaggactctgcgggagttgcgtacctgagcgcggccaccggccgacccatcgcgacaccggcgccctatcgtacctgtgaaggctttccgtgcacatcgggtcaccagtgaaaagttgtgcagagcccaatgagcttcatttccaagctgccacatctctgcctcctgcgtagcatccggaaacatgtgcccctacatcaggaaattcatggcaagggtgagcgcctcgggtgaggagtgctgctgcttgggtcgaagttgccccatcagcctggagggaagggtggagccatga