

KIAA1143-KIAA1143 Gene

Reference: PTXBC008468

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein KIAA1143		
For research use only.			

GENE DESCRIPTION	KIAA1143
GENE ID	57456
ORIGIN SPECIES	Human
ACCESSIONS	BC008468
VECTOR	pDONR223
NCBI SYMBOL	KIAA1143
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcaagcgggaaccagggtatcgtacgtgcggccagccgagccggcgtttctgcccgttcaaggaaagggtcggctacaggagggaacccacgtagagactaagagaattcagccctcagccccagatgaagatggggatcacagtgacaaagaagatgaacagcctcaagt
gggtgttttaaaaaaggagacctgtcagttgaagaagtcagtaaaattaaagcagaataaaggctgccaagcagatgaagaaccaactccagccgatggaagaatcatatcgaaccagtcagcatccctcagatgaaaaatttcaggtttaacagcagctcaaaaaagaag
aagccaaatgaagatgaagtaatacaggactcgggtcaaaaaaactcacaaaaacaattaaaaatagtagctccttctttgacaacgaagatgaaatgagtaa