



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CCDC34-coiled-coil domain containing 34 Gene

Reference: PTXBC008496

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L15; NY-REN-41; RAMA3; coiled-coil domain-containing protein 34; NY-REN-41 antigen; renal carcinoma antigen NY-REN-41; coiled-coil domain containing 34		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 34
GENE ID	91057
ORIGIN SPECIES	Human
ACCESSIONS	BC008496
VECTOR	pDONR223
NCBI SYMBOL	CCDC34
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]