

IMMP2L-IMP2 inner mitochondrial membrane peptidase-like (S. cerevisiae) Gene

Reference: PTXBC008497

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IMMP2L-IT1; IMP2; IMP2-LIKE; mitochondrial inner membrane protease subunit 2; IMP2 inner mitochondrial membrane peptidase-like; IMP2 inner mitochondrial membrane protease-like; inner mitochondrial membrane peptidase 2 like; inner mitochondrial membrane pe		

For research use only.

GENE DESCRIPTION	IMP2 inner mitochondrial membrane peptidase-like (S. cerevisiae)
GENE ID	83943
ORIGIN SPECIES	Human
ACCESSIONS	BC008497
VECTOR	pDONR223
NCBI SYMBOL	IMMP2L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcacagtcacagggtgggtgaaagatacatcaaggcctttgtaaaggctcttctgtggcggtgcctgtggcagtgacttcttggatcgggtcgctgtgtggcaagtagaaggagcatcgcagccttcttgaatcctggggggagccagtcacatctgatgtgtgctttgaa
ccactggaaagtgaggaaatttgaagtacaccgtggtgacattgtatcattgtgtctcttaaaaccagaacagaagatcattaaagagagtgattgtcttgaaggagatattgtcagagatgggagaaaactgaaaaggatatga