

C4orf34-chromosome 4 open reading frame 34 Gene

Reference: PTXBC008502

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C4orf34; small integral membrane protein 14		
For research use only.			

GENE DESCRIPTION	chromosome 4 open reading frame 34
GENE ID	201895
ORIGIN SPECIES	Human
ACCESSIONS	BC008502
VECTOR	pDONR223
NCBI SYMBOL	C4orf34
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagaaggtggattgatccctgtgaatgtgttgcctcatgaacatgcaatgagaagactgatcaatctgttacggcagtcctcactgcacagacacagagtgtcttcagggaattaccgggaccctctggatgaatggcatcagtggtacaatgatcttgtagcctggatggtta
ttgcattgatctgttcttactgagacctcctaataagaggatccagctacctggaaagccaaccagtcctcataatggacaagatccaccagctcctctgttgactaa