

HLA-DQA1-major histocompatibility complex, class II, DQ alpha 1 Gene

Reference: PTXBC008585

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CELIAC1; DQ-A1; HLA-DQA; HLA class II histocompatibility antigen, DQ alpha 1 chain; DC-1 alpha chain; DC-alpha; HLA-DCA; MHC HLA-DQ alpha; MHC class II DQA1; MHC class II HLA-DQ-alpha-1; major histocompatibility complex, class II, DQ alpha 1		

For research use only.

GENE DESCRIPTION	major histocompatibility complex, class II, DQ alpha 1
GENE ID	3117
ORIGIN SPECIES	Human
ACCESSIONS	BC008585
VECTOR	pDONR223
NCBI SYMBOL	HLA-DQA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatcctaaacaaagctctgatctggggcccttgccctgaccacgtgatgagccctgtggaggtgaagacattgtggctgaccacgtcgcccttatgggtgtaaacctgtaccagctctacgggtccctctggccagtacacccatgaattgatggagatgagcagttctacgtggacct
ggggagggaaggagactgtctgtgtttgctgtctcagacaatttagatttgacccgcaattgacatgacaacacatcgctgtcctaaacataacttgaacagctctgattaaacgtcccaactctaccgctgctaccaatgagggtcctgaggtcacagtgttttccaaagtctcccgtagacactg
ggtagcccaacatcctcatctgtctgtggacaacatcttctcctgtgtgtaacatcacatggctgagcaatgggcactcagtcacagaaggtgtttctgagaccagcttctctccaagagtgtatcctcttcaagatcagttacccacccctcctcctctgtgtaggagattatga
ctgcaaggtggagcactggggcctggacaagcctcttctgaaacactgggagcctgagattccagccctatgtcagagctcacagagactgtgtctgcgccctgggattgtctgtggcctctgggcattgtgtgggcactgtttcatcatccgagggcctgcgttcagttgtgtcttc
cagacaccaaggcccttgta