

CCL5-chemokine (C-C motif) ligand 5 Gene

Reference: PTXBC008600

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	D17S136E; SCYA5; SIS-delta; SISd; TCP228; eoCP; C-C motif chemokine 5; T-cell specific protein p288; beta-chemokine RANTES; chemokine (C-C motif) ligand 5; eosinophil chemotactic cytokine; regulated upon activation, normally T-expressed, and presumably se		

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) ligand 5
GENE ID	6352
ORIGIN SPECIES	Human
ACCESSIONS	BC008600
VECTOR	pDONR223
NCBI SYMBOL	CCL5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggtctccgcggcagccctcgtgtcatctcattgctactgccctctgcgtctcgtatctgctcccatattctcggacacacacccctgtcttggctacattgcccgccactgccccgacacatcaaggagtattctacaccagtggcaagtgtccaaccagcagt
cgctctttgtcaccgaaagaaccgccaagtgtgtgccaaccagagaagaatgggttcgggagtacatcaactctttggagatgagctag