

CCL13-chemokine (C-C motif) ligand 13 Gene

Reference: PTXBC008621

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CKb10; MCP-4; NCC-1; NCC1; SCYA13; SCYL1; C-C motif chemokine 13; CK-beta-10; chemokine (C-C motif) ligand 13; monocyte chemoattractant protein 4; monocyte chemotactic protein 4; new CC chemokine 1; small inducible cytokine subfamily A (Cys-Cys), member 1		

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) ligand 13
GENE ID	6357
ORIGIN SPECIES	Human
ACCESSIONS	BC008621
VECTOR	pDONR223
NCBI SYMBOL	CCL13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaagtctctgcagtgcttctgtgctgtgctcatgacagcagcttcaacccccaggagctgctcagccagatgcactcaacgtccatctacttgcctcacatttagcagtaagaagatctcttcagagagctgaagagctatgtgatcaccaccagcaggtgtccccaagg
ctgtcatcttcagaacaaactgggcaaggagatctgtgctgacccaaggagaagtgggtccagaattatataatgaacacctgggccggaagctcacacctgaagactga