

CEP290-centrosomal protein 290kDa Gene

Reference: PTXBC008641

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	3H11Ag; BBS14; CT87; JBTS5; LCA10; MKS4; NPHP6; POC3; SLSN6; rd16; centrosomal protein of 290 kDa; Bardet-Biedl syndrome 14 protein; CTCL tumor antigen se2-2; Meckel syndrome, type 4; POC3 centriolar protein homolog; cancer/testis antigen 87; centrosomal		

For research use only.

GENE DESCRIPTION	centrosomal protein 290kDa
GENE ID	80184
ORIGIN SPECIES	Human
ACCESSIONS	BC008641
VECTOR	pDONR223
NCBI SYMBOL	CEP290
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccatttcaagattgcagctctcaaaaaagttagataatagtttcttctgaactagaactggcctaataacagtcacatgaactgactgctaagtacagggacatcttgcataaaagataatgctgttcaagaacaagtaacttgaacacctggagtgtaaaacatctccta
aaagaacaagtggagtgctataataaagaactggagattaccaaggaaaaacttcacactattgaacaagctgggaacaggaaactaaattagtaattgaatctagcattggaaggaagcaagaatacaataaccaacagtgacattgttccatttcaaaaaataactatgctggaatga
aggaattaaatgaaagcagcgggctgaacattgtcaaaaatgtatgaacacttacggacttcgttaagcaaatggaggacgtaattttgaattggaaaccaaatttctgaggtttga