

MGC9913-hypothetical protein MGC9913 Gene

Reference: PTXBC008651

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein MGC9913
GENE ID	386759
ORIGIN SPECIES	Human
ACCESSIONS	BC008651
VECTOR	pDONR223
NCBI SYMBOL	MGC9913
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatcgaggagacacagactggacctgctacgctctctggagtgatcagactgatgatcaccaaccaactcattccggataaggaagaagagagtgacacttctcagtggtttcaaccctactctgcatcttaagacactgtgtacaacgttgacactgtgcagg
atgatgccacttcatcttggatgctaactgccaatgttgactctga