

GNGT2-guanine nucleotide binding protein (G protein), gamma transducing activity polypeptide 2 Gene
Reference: PTXBC008663

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	G-GAMMA-8; G-GAMMA-C; GNG8; GNG9; GNGT8; guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-T2; G protein cone gamma 8 subunit; G-gamma-9; g gamma-C; gamma-T2 subunit; guanine nucleotide binding protein (G protein), gamma transducing activity		

For research use only.

GENE DESCRIPTION	guanine nucleotide binding protein (G protein), gamma transducing activity polypeptide 2
GENE ID	2793
ORIGIN SPECIES	Human
ACCESSIONS	BC008663
VECTOR	pDONR223
NCBI SYMBOL	GNGT2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccaggatctcagcgagaaggacctgttgaagatggaggtggagcagctgaagaaagaagtgaaaaacacaagaattccgattfccaagcgggaaaggaaatcaaggagtaoigtgaggcccaagcagaaacgatcctttctcaaaggcatccctgaggacaagaatccc
ttcaaggagaaagggtgctgtctgataagctga