

BIRC5-baculoviral IAP repeat-containing 5 Gene

Reference: PTXBC008718

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	API4; EPR-1; baculoviral IAP repeat-containing protein 5; apoptosis inhibitor 4; apoptosis inhibitor survivin; survivin variant 3 alpha; baculoviral IAP repeat containing 5		

For research use only.

GENE DESCRIPTION	baculoviral IAP repeat-containing 5
GENE ID	332
ORIGIN SPECIES	Human
ACCESSIONS	BC008718
VECTOR	pDONR223
NCBI SYMBOL	BIRC5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtgccccgacgttgccccctgctggcagccctttctcaaggaccaccgcattctacattcaagaactggcccttttgagggtgctgcctgcacccggagcggatggcggaggtggttcacactgccccactgagaacgagccagactggccccagtgtttcttctgctt
caaggagctggaaggctggagccagatgacgacccatagaggacataaaaaagcattctcgtggtgctgttcttctgtcaagaagcagtttgagaattaacccttggtgaattttgaaactggacagagaaagccaagaacaaaattgcaaaggaaaccaataagaag
aaagaatttgaggaaactgcgaagaagtgcgctggtccatcgagcagctggctgcatggattga