

TEX264-testis expressed 264 Gene

Reference: PTXBC008742

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZSIG11; testis-expressed sequence 264 protein; testis expressed gene 264; testis expressed 264		

For research use only.

GENE DESCRIPTION	testis expressed 264
GENE ID	51368
ORIGIN SPECIES	Human
ACCESSIONS	BC008742
VECTOR	pDONR223
NCBI SYMBOL	TEX264
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggacctgctactctggcctgattggggcctgactctcttactgctgctgacgtgctgcttggccgggtactcagggtactgctgggtggaagtgaagtgtggtcaccctccatccgcaacgtcactgtggcctacaagtccacatggggtctatggtgagactggg
cggcttttactgagagctgcagcatctctcccaagctccgctccatcgtgtctactatgacaaccccatggtgctccctgataaagtgcctggtggtggcagctcctgagtgaggtaagtcgcctccctgagctcagctctaccagaaattggcttcaagggtg
tctccttccggcaccagccatgtgtgacagcccttccctacaccacattctgtccatctggctggtgctacccgctgtccatctgcttggacacctacatcaaggagcgggaagctgtgtgcttccctcggctggagatctaccaggaagaccagatccattcatgtgccact
ggcacggcaggagacttctatgtgctgagatgaaggagacagagtgaaatggcggggcctgtggaggccattgacacccaggtggtggcagcagagctgacacaatgagtgacacagagttctgtaagcttggagtgagccctggcagcggggagacttcagctgccacact
gtcaccctggggcagcagccgtggctggatgacgggtgacacccgagcagcagctacagcagtcagtgccagcggctctcttggaggagctggacttgaggcgaggggcccttaggggagtcacggctggaccctgggactgagccctggggactaccaagtgg
ctctgggagccactgccccctgagaaggcaaggagtaa