



Reference: PTXBC008743

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ESP-2; HED-2; zyxin-2		
For research use only.			

GENE DESCRIPTION	zyxin
GENE ID	7791
ORIGIN SPECIES	Human
ACCESSIONS	BC008743
VECTOR	pDONR223
NCBI SYMBOL	ZYX
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]