



Reference: PTXBC008747

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CAB1; es64; stAR-related lipid transfer protein 3; MLN 64; START domain-containing protein 3; StAR-related lipid transfer (START) domain containing 3; metastatic lymph node gene 64 protein; metastatic lymph node protein 64; steroidogenic acute regulatory		

For research use only.

GENE DESCRIPTION	StAR-related lipid transfer (START) domain containing 3
GENE ID	10948
ORIGIN SPECIES	Human
ACCESSIONS	BC008747
VECTOR	pDONR223
NCBI SYMBOL	STARD3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]